Flow Visualization Equipment and Facilities 09/10/24

MCEN 4151-5151/ ATLS 4151

Flow Visualization: The Physics and Art of Fluid Flow

Here is a list of flow facilities; equipment for checkout is listed below. Unless otherwise specified, most of this stuff is in my lab, and I will check it out to you. MWF afternoons are the best time to catch me.

ITLL equipment:

The Consumables Closet, ITLL 151 near the Ball Machine (north end of upper level), has all kinds of free stuff to make small fluids apparatuses from. Check with ITLL Launchpoint (south end of upper level) if you don't see what you want.

Idea Forge: The Idea Forge (east end of Fleming) also has a huge assortment of free parts for DIY setups; glassware, plexi, pumps, plumbing, fans etc..

FLOW FACILITIES: AIR

Facility	Lighting	Visualization	Phenomena	Access
Vortex ring	Try projector	Stage fog	Vortex rings,	JH Lab
generators;	for light sheet,		symmetric and	
zeroblaster, or	or strobe		asymmetric.	
timed generator (needs a little work).				
,			Brynne Sutton, Emrys Hall, Thomas King,	
			Bethany Rotherham FV2003	
Laser sheet/fog	Built in rotating	Built-in stage	Turbulent jet	JH Lab
Desk toy	mirror and	fog generator	cross section	
	green laser		and room air	
	pointer		turbulence/mixi	
			ng	

Medical	Strobe for	Dry ice vapor ¹	Jet flows,	JH has
nebulizer,	volume vis	humidifiers,	positive	nebulizers,
ultrasonic		steaming pots,	buoyancy	humidifier, fog
humidifier		medical	convective flow	machines
		nebulizers	<u>Ultrasonic</u>	
		(<\$5) ²	atomizer (same	
		Fog generators	·	

¹ Dry ice is solid carbon dioxide. Do not seal into a container, let it breathe. Handle with extreme care; it can freeze flesh and displace breathable air. Cover with hot water for best effect, otherwise a water ice shell will form.

² Do not nebulize oils (i.e. canola) without use of a proper respirator or aerosol filter mask: oil coated lungs define pneumonia and asphyxiation.

2 small (4 inch diameter) schlieren systems for home checkout. See Guidebook.	Single LED headlight works well.	Schlieren: Light bent by η gradients Could do stereo with 2 small systems	as humidifier) flow Convective flows from warm/hot objects: hands, candles, hair dryers (turbulent jet). You may need time to make your own color stops. Can be used in water too.	Colleen Stroud FV 2004
Reuben's Tube	Flame	Flame length represents pressure.	Standing wave resonance in a pipe, excited by a loudspeaker on the end.	Maybe? These come and go.

FLOW FACILITIES: LIQUIDS

Facility	Lighting	Visualization	Phenomena	Access
Flume	LED panel	Surface flow or food dye	Open channel	ITLL.
		injection. Must change	flow, 8 ft long	Reservati
		water after	x 6inch wide	on and
			x 1 ft deep.	maybe
Tanana I adhlassa Tira			Flow around	short
Tanner Ladtkow, Tim Read FV 2006			obstacles	training
				required.
				Ask at
Melissa Talmage, Nigel				Launch-
Gorbold, Lok Kin lee,				point
Christopher McCray, Taylor				Or
Simonson FV2006				Rachel.sh
				arpe@col
- 11				orado.edu
Small water	Includes	Bubbles Dye, rheoscopic	Designed for	JH lab.
tunnel for	water pump	fluid, paint, anything safe	object wakes	Requires
checkout; 3'	for	for drains		setup.
long, 2'deep	circulation			

Large Fish	Strobe or	Food coloring.	Short jets,	Can set
Tank. 50 gal.	work lights	Be sure to bleach water	vortex rings,	up in
(50 gal)		clean afterwards	boundary	ITLL
			layers	
Several small (10 gal) Fish Tanks, larger fish tank, pumps available too. For checkout.	Strobe, laser sheets	Food coloring, alumina powder, cornstarch particles; anything you are willing to put down your own drain.	Short jets, vortex rings, boundary layers Steady vertical vortex (from stirring machine) Small ring generators	
Hele-Shaw cell	Work light	Food coloring of	available. Saffman-	
Ticic-Silaw cell	or bounced	detergent, corn syrup,	Taylor	
Transparent glass plate Leav viacous thaid Nigh viacous thaid	strobe	water, etc	instability	
Singh, Akhileshwar, Yogesh Singh, and Krishna Murari Pandey. "Viscous Fingering Instabilities in Radial Hele-Shaw Cell: A Review." Materials Today: Proceedings, 10th International Conference of Materials Processing and Characterization, 26 (January 1, 2020): 760–62. https://doi.org/10.1016/j.matpr.2020.01.022.		https://www.flowvis.org/2015/06/01/bl ue-milk-and-then-yellow-milk-were- injected-into-clear-detergent-in-a-hele- shaw-cell-to-create-the-saffman-taylor- fingering-instability/		

Glitter Track	LED or other	Glitter (Pearl-	Wake and wave	JH. Would
6 foot X 3 inch	worklights	Ex), Pearl Swirl	phenomena	benefit from
black PVC half		or pearlescent		small recirc
tubes		shampoo		pump.
Liquid Desk		Built in	Various,	JH office. An
Toys: lava			including low-	assortment of
lamp, vortex			order	dynamic desk
lamp, drip			turbulence,	toys that have
timers, sparkly			wakes, droplet	fluid motion.
fluid in balls,			motion	
etc.				
Blackstock	Has polarized	Streaming	Cylinder wake	Prof. Hertzberg.
Rheoscopic	light setup	birefringence		Also have extra
Fluid cell				fluid available,
				but apparatus

				must be very
				clean; no salts.
Ferrofluid	Normal studio	Move it with	Magnetic field	JH. Bring small
Regular and	lighting	magnets. You	lines show the	container.
mag clutch		need to supply	'normal field	Impossible to
101917		magnets.	instability'	clean up spills.
THE-CLUBAC ORY				Will stain
fluid:			Images with	anything.
			just a magnet	Nontoxic,
Ferrofluid			under a pool of	though.
Climbs			ferrofluid	
http://vimeo.co			showing the normal field	
<u>m/55136676</u>			instability are	
David Oakley, Peter Davis, Kerylyn Lay,			not allowed.	
Jakob Anderegg,			See Team Third	
Brayden Hass. 2012			2023 to see	
			why.	
Ferrofluid Flies				
Up				
http://vimeo.co				
m/55075720				
Brayden Hass, Jakob Anderegg, Peter Davis,				
Kerylyn Lay, David				
Oakley 2012				
Glycerin			Iridescent soap	JH lab. Mix
			bubbles and	with soap
			films.	solutions to
				extend soap
				film life
Droplet Splash	Has dedicated	Reflection and	Worthington	JH
System. <u>Has</u>	Nikon camera	refraction of	jets, crown	3
excellent	and strobes.	fluids	splashes	
documentation	Requires			
	Android phone for synch app			http://www.flow
	Tor synon app			vis.org/2016/09/
				11/worthington-
				<u>jet-of-first-drop-</u> <u>collides-with-</u>
				second-drop/
	<u> </u>	1	1	3000Hu-urup/

Small Equipment Checkout

Equipment	Location	Notes

Stage fog generator (cooled)		JH		Fog is nontoxic water-based glycol solution. Can leave residue, and may trigger smoke alarms in high concentrations.
Stage fog ge	Stage fog generator, (small)			Buy at Lowe's or Target in late September, \$30 Is basically a vape system
		***		but no human required.
City models,	tamous ~ 1 ft ² . 1:5000	JH		Use with water or air. Supply your own Godzillas.
Ultrasonic h		JH		Makes acoustic streaming
				fountain.
Glass sheets,	tempered	JH		28" x 10 to 17 inch. 3
				sheets. Play with surface
				tension effects
		CAMERAS	and LENSES	
ITLL Electro	nics Shop. See J	onah Spicher	<jonah.spicher@< td=""><td>ocolorado.edu></td></jonah.spicher@<>	ocolorado.edu>
Mfg	Тур	е	Model	Link
Sony	Camcorder 4K F	IDR	HXR-NX80	https://pro.sony/ue_US/hxr-
JVC	Camcorder 1920	O HD	GC-PX100BU	https://uk.jvc.com/microsite px100/index.html#SPECIFIC
Kron Technologi es	Slo-Mo 1280x1024 @ 1,069FPS (max res) 40,413 Frames Per Second (FPS) at lower resolution.		Kronos 1.4 High Speed Camera	https://www.krontech.ca/pronos-1-4-high-speed-camera
Panasonic	Mirrorless. 20 Mpx, 4K video		Lumix GH5	https://shop.panasonic.com gh5m2-mirrorless-camera-b
Canon	DSLR. 15 Mpx, HD video, ISO 128,000		EOS Rebel T1I	https://global.canon/en/c- museum/product/dslr801.ht
GoPro	4K video max, wide angle lens, 240 fps at low res		Hero3	https://gopro.com/content/cont
Canon Eos T kit lens, 2 ba	3i DSLR with tteries	JH		
	lapse trigger	JH		

Canon extension tubes (for	JH	
cheap lenses, no electronic		
pass thru)		
Nikon extension tubes	JH	
Nikon 24 mm wide angle	JH	
lens		
Nikon 50 mm lens	JH	
Nikon macro lens 102 mm	JH	Manual only
Closeup Lenses: +1, 2, 4 in	JH	
58 mm dia, +2,+3 in 72 mm		
dia.		
Stereo cameras (film) and	JH	
slide bar		
FLIR Infrared thermal	JH	
imaging camera. For		
iPhones only		
•		
	LIGHTING	
Dimmable Bi-color 660	ITLL has one, JH has two	
LED Video Light		
(continuous) with Barndoor		
and 6.5 feet Light Stand,		
Godox VING V860IIN	JH	Good for remote and
TTL Li-Ion Flash with		multiple triggers
X1T-C TTL Trigger Kit		
for Nikon		
Godox VING V860IIC TTL	JH	Good for remote and
Li-Ion Flash with X1T-C		multiple triggers
TTL Trigger Kit for Canon		
Cameras		
Yongnuo YN560-IV	JH	Can be triggered by flash on
Speedlite Flash with		your camera
Manual and Slave Control"		
Umbrella reflectors	JH	2 on stands, one short/table
		mount. Comes with weak
		CFL lights
24" (60cm) 5-in-1 Disc	JH	Translucent, Silver, Gold,
Light Reflector with Bag -		White and Black
Sunpak Auto 383 Flash	JH	
(strobe) unit & 25' pc cable		
CW 1 watt blue LED laser	JH	Serious safety training
		required
Party strobe	JH	
500 W work lights, several	ITLL, JH	Incandescent. Produces a lot
sets		of heat as well as light.
· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·

Small LED worklight pair	JH	
North Star video lights (2),	Idea Forge	
cooled		
	MISC	
Gretag-Macbeth/X-Rite	See Prof. Hertzberg	For color calibration of
Eye-1 Spectrophotometer		monitors, cameras, printers
		and projectors.
Large black backdrop (8	Idea Forge	Lots of fun to fold back up.
foot square), integral stand		
Small white table-top tent,	Idea Forge	Provides diffuse white light
$\sim 2 \text{ ft}^3$		and control of reflections
black velvet	JH	Small pieces
Assorted tripods	JH	
LP Turntable	JH	For study of rotating flows